

# TRIA™

## Handheld Doppler










- **Sensitec Signal Processing™**  
*SSP™ provides pure, clear sound for greater accuracy in patient diagnosis*
- **Blended Beam™ Technology**  
*Groundbreaking dual-crystal system designed for rapid signal acquisition*
- **Accu-Rate™ Digital Display**  
*NEW! Easy-to-read Accu-Rate™ Digital Heart Rate Display calculates beats per minute to save counting and calculating time*
- **Excellent Value**  
*Tria™ is a durable and elegant Doppler at an affordable price*
- **Flexible**  
*Interchangeable Fetal & Vascular Probes and optional **Waterproof Fetal Model** provide the versatility you need in all clinical settings*

*Tria* **Doppler System**

MEDASONICS//®

# TRIA™ Doppler Ultrasound Stethoscope®

## ORDER INFORMATION

	Description	Order Number	Frequency	Contents - Applications / Specifications
	TRIA Fetal Doppler with Interchangeable Probe or with Dedicated Probe	T334 T345*	2.3 MHz 2.3 MHz	<b>Includes:</b> Doppler probe, integrated speaker unit, gel, battery, soft carrying case. <b>Applications:</b> Detect fetal life early in pregnancy. Assess the rate and rhythm of the fetal heart.
	TRIA Fetal Doppler with Accu-Rate™ Heart Rate Display and Interchangeable Probe or and Dedicated Probe	T334-D T345-D*	2.3 MHz 2.3 MHz	<b>Includes:</b> Doppler probe, integrated speaker/heart rate display unit, gel, battery, soft carrying case. <b>Applications:</b> Detect fetal life early in pregnancy. Assess the rate and rhythm of the fetal heart.
	TRIA Vascular Doppler with Interchangeable Probe or with Dedicated Probe	T335 T346*	5.0 MHz 5.0 MHz	<b>Includes:</b> Doppler probe, integrated speaker unit, gel, battery, soft carrying case. <b>Applications:</b> Locate vessels and assess vessel patency, calculate ankle-brachial index, measure systolic blood pressure.
	TRIA Doppler with Interchangeable Fetal & Vascular Probes w/o heart rate display w/Accu-Rate™ Heart Rate Display	T344 T344-D	2.3 & 5.0 MHz	<b>Includes:</b> Doppler probes (2), integrated speaker unit (optional heart rate display), gel, battery, hard and soft carrying cases. <b>Applications:</b> See above Fetal and Vascular application descriptions.
	TRIA Fetal Doppler with Waterproof Fetal Probe w/o heart rate display w/Accu-Rate™ Heart Rate Display	T336* T336-D*	2.3 MHz	<b>Includes:</b> Waterproof Doppler probe, integrated speaker unit (optional heart rate display), gel, battery, soft carrying case. <b>Applications:</b> Detect fetal life early in pregnancy. Assess the rate and rhythm of the fetal heart.
	Fetal Probe	P338	2.3 MHz	<b>Includes:</b> Fetal Doppler probe, gel.
	Vascular Probe	P339	5.0 MHz	<b>Includes:</b> Vascular Doppler probe, gel.

**Additional Accessories:** Secure your TRIA Doppler with a special "theft-proof" bracket that is available separately or mounted on a heavy duty rollstand. The SH3A Stethoscope headset may also be used with TRIA. Contact your Sales Representative for product information, part numbers and pricing.

**\*Note:** TRIA Dopplers with dedicated probes are available with a longer cord (at no extra charge) for use with our Rollstand or Bracket. Order a TRIA-LC (i.e., T345-LC) when you order a Doppler to go with your Rollstand or Bracket.

## TRIA™ SPECIFICATIONS

<b>Ultrasonic Frequency</b>	Fetal Doppler: 2.3 MHz Waterproof Fetal Doppler: 2.3 MHz Vascular Doppler: 5.0 MHz
<b>Outputs</b>	One 3.5mm jack to connect to a headset or audio recording device
<b>Controls</b>	Volume control and push button ON/OFF button
<b>Power Source</b>	One 9V alkaline battery
<b>Weight</b>	8.75 oz. (245 gm)
<b>Dimensions</b>	5.25 x 2.6 x 1.75 in (13.3 x 6.6 x 4.5 cm)
<b>Construction</b>	Rugged ABS polycarbonate case
<b>Service</b>	CooperSurgical, Inc. maintains a full service facility which can promptly repair products returned to the factory.
<b>Warranty</b>	MedaSonics Dopplers are covered by a 5-year warranty on parts; 1-year warranty on labor. <i>Call Customer Service for details.</i>

# MEDASONICS//®

## COOPERSURGICAL

95 Corporate Drive, Trumbull, CT 06611

800-227-8076 • 203-601-5200

Fax: 800-262-0105

[www.coopersurgical.com](http://www.coopersurgical.com)