

Preemie Pouch®


New Products!

1-piece and 2-piece ostomy pouching system

Reinforced with double hydrocolloid thickness around the stoma

New wider pouch seals offer extra security against leaks

Designed for high output and difficult to pouch infants



Catalog# 7OS2422

2-piece ostomy pouch without plastic rings to snap together. Just gently place the daisy shaped hydrocolloid with double thickness of hydrocolloid around the stoma and simply adhere the pouch. Particularly suitable for infants with constant discharge such as high ileostomies and jejunostomies.

Catalog# 7OS1422

1-piece ostomy pouch with double thickness of hydrocolloid around the stoma to extend the wear time for high output infants. Particularly suitable in ileostomies, jejunostomies and difficult colostomies. Also great for extra security in all neonatal ostomy patients.


Two new unique ostomy pouches specifically developed to meet the need of the neonatal ostomy patient - from the experts in neonatal ostomy care

Incutech Inc.

P.O. Box 1608

Kernersville, NC 27285-1608

Phone: 800-699-4232 Fax: 336-993-4883 e-mail: info@incutech.com



ORDER INFORMATION

Preemie Pouch can be ordered directly from Sun Medical by calling Customer Service at

800-678-6633
or fax
817-640-1840

Cat. Nr.	Product Description
	1-piece and 2-piece neonatal, drainable, clear ostomy pouch with flexible hydrocolloid skin barrier. All pouches can be cut to desired length with a maximum of 12". 10 pouches per box.
7OS1012	1½" round barrier without starter hole
7OS1016	2" round barrier without starter hole
7OS1020	2½" round barrier without starter hole
7OS1222	2¾" daisy shaped barrier without starter hole
7OS1228	3½" daisy shaped barrier without starter hole
7OS1112	1½" round barrier with ¼" starter hole
7OS1116	2" round barrier with ¼" starter hole
7OS1120	2½" round barrier with ¼" starter hole
7OS1322	2¾" daisy shaped barrier with ¼" starter hole
7OS1328	3½" daisy shaped barrier with ¼" starter hole
NEW! 7OS1422	2¾" daisy shaped barrier without starter hole and with double thickness hydrocolloid around the stoma (1 ½" diameter)
NEW! 7OS2422	2¾" daisy shaped barrier without starter hole and with double thickness hydrocolloid around the stoma (1 ½" diameter) ***** 2-piece pouching system *****

Distributed by:

Sun Medical, Inc.
1179 Corporate Drive West, Suite 100
Arlington, TX 76006