



## Kentec Medical's "Accurate" Line of Products

For 20 Years, Kentec Medical has cared for your nursery patients with our "Accurate" Line of disposable products. We are pleased to announce the addition of a complete line of disposable and reusable temperature probes. The Accutemp-Probe was created for the sole purpose of providing optimal temperature maintenance of infants during hospital therapy. Our Accutemp-Probes and Insulated Thermal Reflective Hearts together offer the best combination to attain Accurate temperature readings.



**NEW!**

### Temperature Probes:

#### Accutemp-Probe™ Disposable Skin Temperature Probes are:

- Available in both disposable & reusable styles
- 100% tested for accuracy
- Designed especially for use with Air Shields, Ohmeda, Hill-Rom and Dräger Incubators and Warmers
- Guaranteed O.E.M. compatible
- Conveniently packaged in boxes of 10 and cases of 6 boxes
- Packaged with samples of both hydrogel and adhesive insulated thermal reflective probe covers

### Probe Covers:



#### Accutemp™ Insulated Thermal Reflective Hearts:

You have a choice! Accutemp gentle grip adhesive or individually packaged Accutemp Plus™ hydrogel for sensitive skin and convenient repositioning. The insulating foam layer keeps the ambient air from affecting the true temperature readings while the gold mylar foil deflects radiant heat. These factors assure that the probe is being warmed by the infant, not the incubator or warmer.

#### Ordering Information:



- # 120 Accutemp Insulated Thermal Reflective Hearts  
Adhesive Box of 120
- # 120P Accutemp Insulated Thermal Reflective Hearts  
Adhesive Box of 120 (Premie Size)
- # 175\* Accutemp-Plus Insulated Thermal Reflective Hearts  
Hydrogel Box of 75
- # 175-6\* Accutemp-Plus Insulated Thermal Reflective Hearts  
Hydrogel Ctn./6 Bxs.

\* Individually packaged hydrogel





# Ordering Information



For use with Air Shields C-100, C-450, C-550 Incubators and PM-78 Infant Warmers.

Catalog #210 Box of 10 Disposable Probes

Catalog #211 Case of 6 boxes of 10 disp. probes.



Catalog #215 Reusable Probe, one each.

For use with Air Shields C-86, Incubator and SC-78 Infant Warmers.

Catalog #220 Box of 10 Disposable Probes

Catalog #221 Case of 6 boxes of 10 disp. probes.

Catalog #225 Reusable Probe, one each.



For use with Ohmeda I.C., CG and C.P. Incubators and I.W.S. Infant Warmers.

Catalog #230 Box of 10 Disposable Probes

Catalog #231 Case of 6 boxes of 10 disp. probes.


Catalog #235 Reusable Probe, one each.

For use with Ohmeda 190 ASC Incubator and NICC Warmer.

Catalog #237 Box of 10 Disposable Probes

Catalog #238 Case of 6 boxes of 10 disp. probes.

Catalog #236 Reusable Probe, one each.



For use with Hill Rom Stabilet<sup>®</sup> 200, 2000 A/B and 300, 3000 A/B Series Neonatal Care Center.


Catalog #240 Box of 10 Disposable Probes

Catalog #241 Case of 6 boxes of 10 disp. probes.

For use with Hill Rom Stabilet<sup>®</sup> CC, 2000C, 1500CC and 1250 Neonatal Care Center.

Catalog #243 Box of 10 Disposable Probes

Catalog #244 Case of 6 boxes of 10 disp. probes.



For use with Dräger 8000 Series Incubator.

Catalog #246 Box of 10 Disposable Probes

Catalog #247 Case of 6 boxes of 10 disp. probes.

Distributed By:



17871 Fitch,  
Irvine, CA 92614-6071  
800-825-5996