

**New**

# LIGHT-SAFE™

BC 565

## Amber Tubing Extension Set

- Greatly reduces hydroperoxide generation in neonatal patients receiving lipids and TPN

- Proven to greatly reduce photo-degradation of light sensitive infusates including Amphotericin B, Adriamycin, Furosemide, Dacarbazine, Nitroprusside, Vitamin A

- Priming volume 2.0 ml

- 60" minibore amber tubing extension set

- Totally Inert PVC free fluid path

For more information on LIGHT-SAFE™ call CODAN US Corporation 1-800-33CODAN

